

仅 供 科 研 使 用

HEPES 液

(1M,pH7.2~7.5)

■ 品牌: Oumarsi/欧玛仕 ■ 规格: 100ml/500ml

■ 货号: OMDCM016-100ml ■ 用途: 仅供科研使用

■ 储存: 2-8°C ■ 保质期: 12 个月



产品

HEPES 在 PH 7.2~7.5 时,具有较好的缓冲能力,终浓度一般为 10~50mmol/L。根据具体实验,一般工作浓度为 10~50mmol/L。

对细胞无毒性,优点是在开放式培养或细胞观察时能维持较恒定的 P H 值。高浓度 HEPES 对某些细胞可能有一定毒性,操作注意避免溶液被微生物污染

注意事项:

2-8℃, 12 个月

相关产品:

货号	产品名称	规格	保存条件
OMDCFM-032	组织保存液	50ml/100ml	-20℃,避 光,24 个 月。
OMDCM-004	高糖 DMEM 培养基	500ml	12 个月
OMDCM-02	D-Hanks 平衡盐溶液(1 ×,无酚红)	500ml	12 个月
OMDPBS-011	Dulbecco's 磷酸盐缓冲 液(DPBS),含钙、镁, 不含酚红	500ml	室温,12 个月。
OMDPBS-02	Dulbecco's 磷酸盐缓冲 液(DPBS),不含钙、 镁,不含酚红	500ml	室温,12 个月。
OMDCM-047	α-MEM 培养基	500ml	2-8°C,12 个月。



货号	产品名称	规格	保存条件
OMDCM-040	McCoy's 5A 培养基	500ml	2-8℃,12
OMDCM-031	DMEM 培养基 (高糖,无酚红,含胎牛血清,含双抗)	500ml	2-8℃,12 个月。
OMDCM038-500ml	DMEM 培养基(无糖,含酚 红)	500ml	2-8℃,12
OMDCM017-100ml	硫酸庆大霉素 (50mg/ml)	100ml	2-8°C,12 个月。
OMDCM016-100ml	HEPES 液 (1M,pH 7.2~ 7.5)	100ml/500ml	2-8°C, 12 个月。
OMDPBS-18	1×PBS(0.0067M,pH 7.0-7.2)	500ml	2-8℃,12
OMDPBS-15	1×PBS(0.01M,pH9.0)	500ml	2-8°C, 12 个月。
OMDPBS-06	20×PBS(0.2M,pH 7.2-7.4)	500ml	2-8°C,12 个月。
OMDPBS-010	1×PBS(0.01M,pH6.8)	500ml	2-8°C,12 个月。
OMDPBS-12	1×PBS(0.01M,pH 8.9)	500ml	2-8°C,12 个月。
OMDPBS-20	1×PBS(0.01M,pH 7.2-7.4)	500ml×10	2-8°C,12 个月。



货号	产品名称	规格	保存条件
OMDCM027-500ml	DMEM 培养基(无糖,含酚	100ml	2-8°C, 12
OMDCM027-300IIII	红,含双抗)	1001111	个月。
OMDCM025-500ml	DMEM 培养基(无糖,无酚	100ml/500	2-8°C, 12
OMDCM023-300IIII	红)	ml	个月。
OMDCM024-500ml	DMEM培养基(无糖,无酚	500ml	2-8°C, 12
OMDCM024 300III	红,含双抗)	3001111	个月。
OMDCM022-500ml	DMEM 培养基(高糖,无丙	500ml	2-8°C, 12
OMDEMOZZ SOOM	酮酸钠,含酚红)	3001111	个月。
OMDCM021-500ml	DMEM 培养基(高糖,无丙	500ml	2-8°C, 12
OMDOMOZI SOOM	酮酸钠,含酚红,含双抗)	3001111	个月。
OMDCM019-500ml	DMEM 培养基(高糖,无酚	500ml	2-8°C, 12
	红)	300 1111	个月。
OMDCM018-500ml	DMEM 培养基(高糖,无酚	500ml	2-8°C, 12
	红,含双抗)	Jouint	个月。
OMDCM017-500ml	DMEM 培养基 (低糖,无酚	500ml	2-8°C, 12
	红,含胎牛血清,含双抗)		个月。
OMDCM016-500ml	DMEM 培养基(低糖,无酚	500ml	2-8°C, 12
	红)	3001111	个月。
OMDCM015-500ml	DMEM 培养基 (低糖,无酚	500ml	2-8°C, 12
	红,含双抗)		个月。
OMDCM014-500ml	RPMI 1640 培养基(无酚	500ml	2-8°C, 12
	红,含胎牛血清,双抗)	0001111	个月。



货号	产品名称	规格	保存条件
OMDCM013-500ml	RPMI 1640 培养基(无 酚红)	500ml	2-8°C,12 个月。
OMDCM012-500ml	RPMI 1640 培养基(无酚 红,含双抗)	500ml	2-8℃,12 个月。
OMDCM-017	Myco-Out 支原体清除 试剂(100X)	10ml	2-8°C,12 个月。
OMDCM-015	Hanks 平衡盐溶液(1×, 无酚红) (pH7.2-7.4)	500ml	室温,12 个
OMDCM017-500ml	DMEM 培养基 (低糖,无酚红,含胎牛血清,含双抗)	500ml	2-8℃,12
OMDCM-030	青霉素-链霉素-两性霉素 B 溶液(100×三抗)	100ml	-20°C,12 个月。
OMDCE-050	青霉素-链霉素溶液(100 ×)	100ml	-20°C,12 个月。
OMDCE-060	胰酶细胞消化液(0.25% 胰酶,含酚红,不含 EDTA)	100ml	-20℃,12
OMDCM-059	胰酶-EDTA 消化液 (0.25%胰酶,含酚红)	100ml/500ml	-20°C,12 个月。
OMDCE-004	胰酶细胞消化液(0.25% 胰酶,不含 EDTA 和酚 红)	100ml	-20℃,12 个月。



货号	产品名称	规格	保存条件
OMDCM-034	胰酶-EDTA消化液(0.25%	100ml/500	-20°C,12 个
	胰酶,不含酚红)	ml	月。
OMDCM-020	胰酶-EDTA消化液(0.05%	100ml	-20°C,12 个
	胰酶,含酚红)		月。
OMDCM-011	胰酶-EDTA消化液(0.05%	100ml/500	-20°C,12个
	胰酶,不含酚红)	ml	月。
OMDCM007-1ml	DAPI(1mg/ml)	1ml	-20°C,12 个
	J,(2g,,		月,避光。
	DMEM (Low Glucose,		2-8°C,12 个
OMDCM-051	with FBS,	500mL	月。
	Penicillin-Streptomycin)		
OMDCM-052	Medium 199 (with FBS,	500mL	2-8°C,12 个
	Penicillin-Streptomycin)		月。
	DMEM (Low Glucose,		2-8℃,12 个
OMDCM-053	with	500mL	月。
	Penicillin-Streptomycin)		
OMD81714-500m	DMEM (Low Glucose)	500mL	2-8°C,12 个
<u>l</u>	•		月。
OMD81711-500m	Medium 199	500mL	2-8°C,12个
l			月。
OMDCM-019	Medium 199 (with	500mL	2-8℃,12个
0 m 0 m 0 m	Penicillin-Streptomycin)		月。



货号	产品名称	规格	保存条件
OMD81717-500ml	RPMI 1640 Medium (with Penicillin-Streptomycin)	500mL	2-8℃,12 个月。
OMDCM-039	Ham's F-12K Medium (with FBS, Penicillin-Streptomycin)	500mL	2-8℃,12 个月。
OMD81708-500ml	Ham's F-12K Medium	500mL	2-8°C,12 个月。
OMDCM-038	RPMI 1640 Medium (with FBS, Penicillin-Streptomycin)	500mL	2-8℃,12 个月。
OMDCM-037	Ham's F-12K Medium (with Penicillin-Streptomycin)	500mL	2-8°C,12 个月。
OMDCM-037	Ham's F-12K Medium (with Penicillin-Streptomycin)	500mL	2-8℃,12 个月。
OMDCM-036	DMEM (High Glucose, with FBS, Penicillin-Streptomycin)	500mL	2-8℃,12 个月。
OMDCM-009	MEM (with Penicillin-Streptomycin)	500mL	2-8℃,12 个月。



货号	产品名称	规格	保存条件
OMDCM-014	DMEM (High Glucose, with Penicillin-Streptomycin)	500mL	2-8°C,12 个月。
OMDCM-026	MEM	500mL	2-8°C,12 个月。
OMD81701-500ml	DMEM/F12 (with FBS, Penicillin-Streptomycin)	500mL	2-8°C,12 个月。
OMDCM-028	MEM (with FBS, Penicillin-Streptomycin)	500mL	2-8℃,12 个月。
OMDPBS-01	10×PBS(0.1M,pH 7.2-7.4)	500mL	2-8℃,12 个月。
OMD81702-500ml	DMEM/F12	500mL	2-8℃,12 个月。
OMDCM-033	DMEM/F12 (with Penicillin-Streptomycin)	500mL	2-8°C,12 个月。
OMDCE-008-500ml	超纯水(无菌)	500mL	室温,12 个 月。
OMDCE-009	碳酸氢钠溶液(7.5%)	500mL	-20℃,12 个月。
OMDCM-013	D-Hanks 平衡盐溶液(1×,含酚红)	500mL	2-8℃,12 个月。